

9853 Pacific Heights Blvd. Suite D. San Diego, CA 92121, USA Tel: 858-263-4982

Email: info@abeomics.com

Monoclonal Antibody to Cytokeratin, pan (Epithelial Marker)(Clone : PAN-CK (Cocktail)

Clone Name : Monoclonal
Clone Name : PAN-CK (Cocktail)
Application : IHC, FACS, WB, IF

Reactivity: Rat
Gene: KRT77
Gene ID: 374454
Uniprot ID: Q7Z794
Format: Purified
Alternative Name: KRT77||KRT1B

Mouse IgG, kappa

Description

Isotype:

Twenty human keratins are resolved with two-dimensional gel electrophoresis into acidic (pI 6.0) subfamilies. This antibody cocktail recognizes acidic (Type I or LMW) and basic (Type II or HMW) cytokeratins, with 67kDa (CK1); 64kDa (CK3); 59kDa (CK4); 58kDa (CK5); 56kDa (CK6); 55kDa (CK7); 52kDa (CK8); 56.5kDa (CK10); 53kDa (CK13); 50kDa (CK14); 50kDa (CK15); 48kDa (CK16); 46kDa (CK17); 45kDa (CK18) and 40kDa (CK19). Many studies have shown the usefulness of keratins as markers in cancer research and tumor diagnosis. KRT-PAN is a broad spectrum anti pancytokeratin antibody cocktail, which differentiates epithelial tumors from non-epithelial tumors e.g. squamous vs. adenocarcinoma of the lung, liver carcinoma, breast cancer, and esophageal cancer. It is useful in characterizing the source of various neoplasms and to study the distribution of cytokeratin containing cells in epithelia during normal development and during the development of epithelial neoplasms. This antibody stains cytokeratins present in normal and abnormal human tissues and shows high sensitivity in the recognition of epithelial cells and carcinomas.

Product Info

purification : Affinity Chromatography

Content: 100 µg in 500 µl PBS containing 0.05% BSA and 0.05% sodium azide. Sodium

azide is highly toxic.

Storage condition: Store the antibody at 4°C; stable for 6 months. For long-term storage; store at

-20°C. Avoid repeated freeze and thaw cycles.

Application Note

Flow Cytometry (0.5-2ug/million cells); Immunofluorescence (1-2ug/ml); Western Blot (0.5-2ug/ml for 2 hours at RT); Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT), (Staining of formalin-fixed tissues requires heating tissue sections i



9853 Pacific Heights Blvd. Suite D. San Diego, CA 92121, USA Tel: 858-263-4982

Email: info@abeomics.com

